

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Short-sleeved Patient Examination Gown

File no.: SPEC-IS8PG
Version-004

ISOL8 HEALTHCARE LTD

Unit 1 Westside, Monavalley Business
Park, Tralee, Co. Kerry, V92 K258,
Ireland

Contact Details

info@isol8health.com
+353 (0) 66 71 21659
www.isol8health.com

Comfortable and high level of dignity

Our single-use short-sleeved patient examination gown is a soft breathable gown which is opaque in colour to offer the wearer a high level of dignity. It closes fully at the rear by way of 2 ties and allows easy examination of the patient by healthcare professionals during patient encounters and surgical or diagnostic procedures. The SMS fabric offers protection against light splashes and mechanical movements. This gown comes in size large which is built to fit most body shapes and sizes and also comes in size XXL for large framed persons.



Product Photos



Product Specification

Product Name	Short-sleeved Patient Examination Gown
Manufacturer Product Codes (MPN)	IS8PG01 (L), IS8PG02 (XXL)
Manufacturer	ISOL8 Healthcare Ltd.
Gown type	Opaque patient examination gown – Single use
Sizes (see size guide below for dimensions)	L (Suitable for most sizes) XXL (Large framed persons)
Fabric	SMS 40gsm (100% polypropylene)
Gown cut	Short sleeved
Neck closure	Ties
Colour	Navy Blue
Sterility	Non-sterile
Ties	Rear ties
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	3 years
Product packaging (inner pack)	10 gowns per (labelled) polythene bag
Product packaging (outer carton)	10 x packs per carton (100 gowns) Outer carton: 600x390x450 mm
Recycling Information	Outer carton is 100% recyclable (cardboard), Inner pack is 100% recyclable (PE), Gowns are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485

Size Guide

Size	Length (cm)	Hem circumference (cm)
L	105	150
XXL	105	165

Product Codes

Product code	Size	Quantity/carton	GTIN (gown)	GTIN (carton)
IS8SG01	L	100	05391541270515	05391541270904
IS8SG02	XXL	100	05391541270522	05391541270911



ISOL8 patient gowns are CE certified class I medical devices intended for hospital patients. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2021 ISOL8 HEALTHCARE LTD
Unit 1 Westside, Monavalley Business
Park, Tralee, Co. Kerry, Ireland, V92 K258
Company No. 687076

info@isol8health.com
+353 (0) 66 71 21659
www.isol8health.com

ISOL/8[®]

*Designed by Healthcare Professionals,
For Healthcare Professionals*