

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



FLUIDSAFE Procedure Mask with Visor

File no.: SPEC-IS8RM03P
Version-003

ISOL8 HEALTHCARE LTD

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35
Ireland

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

Fluid resistant and extra soft

Our **FLUIDSAFE** anti-fog procedure mask is manufactured from high quality fabric with a super soft silk-like inner material providing excellent comfort to the wearer. The anti-fog foam strip prevents fogging and the adjustable face visor provides enhanced protection to the upper face against fluids and particulates. Our type IIR masks provide the highest barrier performance and their fluid resistant properties offer protection against splashing during procedures.

Our procedure masks are most suited for use in clinical environments for patient encounters, examinations and minor procedures. Elastic earloops are easy to apply and provide a comfortable fit.



FLUIDSAFE™

Packaging Information

Product packaging (inner pack)	50 masks per dispenser box Dispenser box: 260x245x80 mm
Product packaging (outer carton)	8 inner boxes per carton (400 masks), Outer carton: 510x280x340 mm, Gross weight: 5.7Kg
Recycling Information	Inner box and outer carton are 100% recyclable (cardboard), Mask is non-recyclable

Product Specification

Product Name	FLUIDSAFE Procedure Mask with Visor
Manufacturer Product Codes (MPN)	IS8RM03P
Manufacturer	ISOL8 Healthcare Ltd.
Mask type	Fluid resistant anti-fog procedure mask with face visor. Suitable for use with spectacles and eyewear.
Compliance with standard	EN14683:2019 Type IIR
Mask size (Length x Width)	175x95 mm
Mask fit	Earloops
Colour	Orange stripes with black UV strip (Other colours available on request)
Sterility	Non-sterile
Mask fabric (3-layers)	Outer: 35gsm Spunbond Polypropylene, Middle: 30gsm Melt-blown Polypropylene, Inner: 20gsm Extra-soft Spunbond
Earloops	Spandex polyurethane coated with nylon thread 160±5mm x 3mm
Nose bar material	Full length aluminium nose piece
Visor/Face Shield	Double-sided anti-fog treated PET (0.1mm thick) 295x125 mm
Anti-fog nose strip	Full length soft anti-fog foam strip
Black UV Strip	40gsm Spunbond Polypropylene
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541270324
GTIN code outer carton	05391541270850



ISOL8 masks are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2023 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**