

# ISOL/8<sup>®</sup>

## SPECIFICATION SHEET

*Designed by Healthcare Professionals,  
For Healthcare Professionals*



**Basics Operating Cap (tie-back)**

File no.: SPEC-IS8HB02

### ISOL8 HEALTHCARE LTD

Unit 4 Joyce House,  
Barrack Square, Ballincollig, Co Cork,  
P31 HW35,  
Ireland

### CONTACT

info@isol8health.com  
IE: +353 (0) 66 71 21659  
GB: +44 (0) 191 816 1966  
www.isol8health.com


# An affordable solution for your operating department

Our range of Basics headwear are all made from a soft, breathable spunbond polypropylene fabric which offers a great combination of both comfort and affordability. Our Basics range has been designed to meet the basic needs of your operating department while offering great value for your money.

This single-use basics operating cap ties comfortably behind the head and offers extended depth allowing greater coverage of the side of the head and is suitable for use by persons with both short and longer hair.



Basics



Basics

## Dimensions

Crown (Length x Width)	26 x 15 cm $\pm$ 1cm
Depth	18cm $\pm$ 1cm
Length/circumference (including ties)	100cm $\pm$ 5cm

# Product Specification

<b>Product Name</b>	Basics Operating Cap (tie-back)
<b>Manufacturer Product Codes (MPN)</b>	IS8HB02
<b>NPC</b>	BWF85012
<b>Manufacturer</b>	ISOL8 Healthcare Ltd.
<b>Product type</b>	Single use operating cap – Extended depth
<b>Sizes (see dimensions below)</b>	One size fits all
<b>Fabric</b>	Spunbond Polypropylene 25gsm
<b>Colour</b>	Blue
<b>Sterility</b>	Non-sterile
<b>Ties</b>	Rear head ties
<b>Latex status</b>	Not made with natural rubber latex
<b>Metal status</b>	Nickel and cadmium free
<b>Shelf-life</b>	3 years
<b>Product packaging (inner box)</b>	100 hats per dispenser box 330x140x185 mm
<b>Product packaging (outer carton)</b>	6 boxes per carton (600 hats)
<b>Recycling Information</b>	Outer carton and inner box are 100% recyclable (cardboard) Hats are non-recyclable
<b>CE certification</b>	MDR 2017/745 – Class I medical device
<b>Manufacturer QMS</b>	ISO 13485
<b>GTIN code inner box</b>	05391541270416
<b>GTIN code outer carton</b>	05391541270980





ISOL8 theatre headwear products are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2023 ISOL8 HEALTHCARE LTD  
Unit 4 Joyce House, Barrack Square, Ballincollig,  
Co Cork, P31 HW35. Company No. 687076

[info@isol8health.com](mailto:info@isol8health.com)  
IE: +353 (0) 66 71 21659  
GB: +44 (0) 191 816 1966  
[www.isol8health.com](http://www.isol8health.com)

# ISOL/8<sup>®</sup>

**Designed by Healthcare Professionals,  
For Healthcare Professionals**

Date of issue: 22/03/23, Version 002