

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Elemental Poly-reinforced Surgical Gown

File no.: SPEC-IS8ER

ISOL8 HEALTHCARE LTD


Unit 4 Joyce House,
Barrack Square, Ballincollig, Co Cork,
P31 HW35,
Ireland

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

High level of performance without compromising your comfort

Our Elemental Poly-reinforced gown (High Performance) provides an affordable solution for procedures where moderate to high levels of fluid are expected. The added reinforcement in the critical zones of the gown offers enhanced protection against strike-through from fluid or blood. Constructed from 43gsm SMMS fabric, it provides high level barrier properties while also offering a high level of breathability to ensure optimal comfort for the wearer. The ergonomic fit with raglan sleeves provides the wearer with added comfort, a greater range of movement and allows for easier donning.



Elemental™

Product Photos



Product Specification

Product Name	Elemental Poly-reinforced Surgical Gown
Manufacturer Product Codes (MPN)	IS8E08 (M), IS8E09 (L), IS8E10 (XL), IS8E11 (XL-Long), IS8E12 (2XL), IS8E13 (2XL-Long), IS8E14 (3XL)
Manufacturer	ISOL8 Healthcare Ltd.
Gown type	Reinforced high-performance sterile surgical gown – Single use only
Compliance with standard	EN13795 - High Performance
Sizes (see size chart below for dimensions)	M, L, XL, XL-Long, 2XL, 2XL-Long, 3XL
Fabric	SMMS 43gsm (100% polypropylene), 35gsm Polyethylene and soft PP breathable reinforcement in critical zones
Gown cut	Raglan sleeves for better unrestricted movement and easier donning
Sleeve Seam	Ultrasonic welded
Neck closure	Hook and loop
Cuffs	Soft knitted cuffs
Anti-static	Yes
Alcohol repellency	Yes
Blood repellency	Yes
Colour	Blue
Sterility	Sterile
Pack contents	2 x sterile scrim towels – 50x40 cm 1 x SMS sterile wrap – 60x60 cm
Ties	Double bonded ultrasonic welded SMS ties, Inner ties: 40±2 cm (both), Outer ties: 75±5 cm and 40±5 cm
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	3 years
CE certification	MDD 93/42/EEC – Class Is medical device Notified body ICIM (CE 0425)
Manufacturer QMS	ISO 13485

Packaging Information

Product packaging (inner pack)	Gown pack (sterile pouch) 340x230x50±20 mm
Product packaging (outer carton)	40 gown packs per carton 600x390x420 mm
Recycling Information	Outer carton is 100% recyclable (cardboard) Inner pack and product are non-recyclable

Size Guide

Size	Length (cm)	Hem circumference (cm)
M	115	145
L	127	150
XL	135	155
XL-Long	150	155
2XL	145	160
2XL-Long	160	160
3XL	150	165

Product Codes

MPN	NPC	Elemental gown type	Qty/ carton	GTIN (gown)	GTIN (carton)
IS8E08	BWK85313	Poly-Reinforced M	40	05391541270089	05391541270607
IS8E09	BWK85316	Poly-Reinforced L	40	05391541270096	05391541270614
IS8E10	BWK85328	Poly-Reinforced XL	40	05391541270102	05391541270621
IS8E11	BWK85329	Poly-Reinforced XL-Long	40	05391541270119	05391541270638
IS8E12	BWK85331	Poly-Reinforced 2XL	40	05391541270126	05391541270645
IS8E13	BWK85324	Poly-Reinforced 2XL-Long	40	05391541270133	05391541270652
IS8E14	BWK85322	Poly-Reinforced 3XL	40	05391541270140	05391541270669



ISOL8 gowns are CE certified class I sterile medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2024 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 17/04/24, Version 007