

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Comfort Surgeon's Hood

File no.: SPEC-IS8HC05

ISOL8 HEALTHCARE LTD
5 Catkin Way, Tindale Crescent
Bishop Auckland, Co Durham
DL14 9TF
United Kingdom

CONTACT
info@isol8health.com
GB: +44 (0) 191 816 1966
IE: +353 (0) 66 71 21659
www.isol8health.com


Extra soft, breathable and sustainable

Our range of Comfort headwear are all made from a lightweight, extra soft 100% viscose fabric made from eucalyptus and birch-pulp, which has an extra soft feel. The fabric is renewable, biodegradable, and sustainable. Air vents in the fabric provide excellent breathability, which along with its superior comfort makes these hats perfect for prolonged use in the operating theatre.

This single-use comfort surgeon's hood offers maximum coverage of the head and neck leaving only the face visible. It is perfect for long procedures or use by persons with long sideburns or long hair.



Comfort



Comfort

Dimensions

Crown (Length x Width)	22 x 21 cm \pm 1cm
Depth (front)	5cm \pm 1cm
Depth (rear)	25cm \pm 1cm
Length/circumference (including ties)	115cm \pm 5cm

Product Specification

Product Name	Comfort Surgeon's Hood
Manufacturer Product Codes (MPN)	IS8HC05
NPC	BWF85014
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Single use surgeon's hood providing maximum coverage to the wearer's head
Sizes (see dimensions below)	One size fits all
Fabric	100% viscose made from eucalyptus and birch-pulp
Colour	Light blue
Sterility	Non-sterile
Ties	Neck ties
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	3 years
Product packaging (inner box)	100 hats per dispenser box 330x140x185 mm
Product packaging (outer carton)	6 boxes per carton (600 hats)
Recycling Information	Outer carton and inner box are 100% recyclable (cardboard) Hats are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541270393
GTIN code outer carton	05391541270966



ISOL8 theatre headwear products are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD
5 Catkin Way, Tindale Crescent, Bishop Auckland,
Co Durham, DL14 9TF, United Kingdom

info@isol8health.com
GB: +44 (0) 191 816 1966
IE: +353 (0) 66 71 21659
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 19/03/25, Version 003