

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Basics Shoe Cover

File no.: SPEC-IS8SC01

ISOL8 HEALTHCARE LTD

Head office:

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland

UK office:

Unit 5 Catkin Way, St Helen
Auckland, Bishop Auckland,
Co Durham, DL14 9TF,
United Kingdom

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

Minimize the risk of contamination

Our shoe covers are intended for use in hospital settings to prevent the transfer of microorganisms and dirt into the clean or sterile environment.

Our Basics shoe cover is ideal for use in minimally contaminated areas where the risks associated with contamination would be low and there is expected to be minimal mechanical wear.

Product Photos



Basics



Basics

Product Specification

Product Name	Basics Shoe Cover
Manufacturer Product Codes (MPN)	IS8SC01
NPC	BWP85007
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Shoe cover
Sizes	One size fits all
Fabric	Spunbond polypropylene 35gsm
Sole	Plain
Fastening	Elasticated rim which fits over the shoe
Sterility	Non-sterile
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner pack)	100 pieces per polythene bag
Product packaging (outer carton)	10 x 100 packs per carton (1000 shoe covers)
Recycling Information	Outer carton is 100% recyclable (cardboard), Inner pack is 100% recyclable (PE), Shoe covers are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541271079
GTIN code outer carton	05391541271123



ISOL8 shoe covers are CE certified class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
GB: +44 (0) 191 816 1966
IE: +353 (0) 66 71 21659
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 10/04/25, Version 004