

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Impervious Shoe Cover (PE)

File no.: SPEC-IS8SC03

ISOL8 HEALTHCARE LTD

Head office:

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland

UK office:

Unit 5 Catkin Way, St Helen
Auckland, Bishop Auckland,
Co Durham, DL14 9TF,
United Kingdom

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

Minimize the risk of contamination

Our shoe covers are intended for use in hospital settings to prevent the transfer of microorganisms and dirt into the clean or sterile environment.

Our PE impervious shoe cover offers the best protection for situations with high contamination risks. The polyethylene material is effective in preventing the transfer of fluids and is suitable if there is expected to be higher mechanical wear.

Product Photos



Product Specification

Product Name	Impervious Shoe Cover (PE)
Manufacturer Product Codes (MPN)	IS8SC03
NPC	BWP85008
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Shoe cover
Sizes	One size fits all
Fabric	30 microns PE (100% Polyethylene)
Sole	Plain
Fastening	Elasticated rim which fits over the shoe
Sterility	Non-sterile
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner pack)	100 pieces per polythene bag
Product packaging (outer carton)	20 x 100 packs per carton (2000 shoe covers)
Recycling Information	Outer carton is 100% recyclable (cardboard), Inner pack is 100% recyclable (PE), Shoe covers are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541271093
GTIN code outer carton	05391541271147



ISOL8 shoe covers are CE certified class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 10/04/25, Version 004