

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Basics Theatre Cap (tie-back)

File no.: SPEC-IS8HB01

ISOL8 HEALTHCARE LTD

Head office:

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland

UK office:

Unit 5 Catkin Way, St Helen
Auckland, Bishop Auckland,
Co Durham, DL14 9TF,
United Kingdom

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

An affordable solution for your operating department


Our range of Basics headwear are all made from a soft, breathable spunbond polypropylene fabric which offers a great combination of both comfort and affordability. Our Basics range has been designed to meet the basic needs of your operating department while offering great value for your money.

This single-use Basics theatre cap ties comfortably behind the head and is most suitable for persons with short hair.


Available in the following colours:

Blue, White, Green, Yellow

Other colours available upon request.



Basics



Basics

Dimensions

Crown (Length x Width)	26 x 15 cm \pm 1cm
Depth	13cm \pm 1cm
Length/circumference (including ties)	100cm \pm 5cm

Product Specification

Product Name	Basics Theatre Cap (tie-back)
Manufacturer Product Codes (MPN)	IS8HB01, IS8HB01W, IS8HB01G, IS8HB01Y
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Single use theatre cap – standard depth
Sizes (see dimensions below)	One size fits all
Fabric	Spunbond Polypropylene 25gsm
Sterility	Non-sterile
Ties	Rear head ties
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner box)	100 hats per dispenser box 330x140x185 mm
Product packaging (outer carton)	6 boxes per carton (600 hats)
Recycling Information	Outer carton and inner box are 100% recyclable (cardboard) Hats are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485

Product variations

MPN	Colour	GTIN (box)	GTIN (carton)
IS8HB01	Blue	05391541270409	05391541270973
IS8HB01W	White	05391541272557	05391541272564
IS8HB01G	Green	05391541272571	05391541272588
IS8HB01Y	Yellow	05391541272588	05391541272601



ISOL8 theatre headwear products are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 16/05/25, Version 004