

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Basics Bouffant Cap

File no.: SPEC-IS8HB06

ISOL8 HEALTHCARE LTD

Head office:

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland

UK office:

Unit 5 Catkin Way, St Helen
Auckland, Bishop Auckland,
Co Durham, DL14 9TF,
United Kingdom

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

An affordable solution for your operating department

Our range of Basics headwear are all made from a soft, breathable spunbond polypropylene fabric which offers a great combination of both comfort and affordability. Our Basics range has been designed to meet the basic needs of your operating department while offering great value for your money.

This single use basics bouffant cap offers a light and spacious design with a round shape allowing full coverage of the wearer's hair. This cap has elastic around the entire circumference which ensures it sits comfortably and securely on the wearer's head. It is suitable for use by person's with both long and short hair.

Available in the following colours:

Blue, White, Green, Yellow

Other colours available upon request.



Basics



Basics

Product Specification

Product Name	Basics Bouffant Cap
Manufacturer Product Codes (MPN)	IS8HB06, IS8HB06W, IS8HB06G, IS8HB06Y
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Single use theatre bouffant cap with an elasticated circumference
Sizes	One size fits all
Fabric	Spunbond Polypropylene 25gsm
Sterility	Non-sterile
Ties	Elasticated circumference
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner box)	100 hats per dispenser box 240x240x220 mm
Product packaging (outer carton)	6 boxes per carton (600 hats)
Recycling Information	Outer carton and inner box are 100% recyclable (cardboard) Hats are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485

Product variations

MPN	Colour	GTIN (box)	GTIN (carton)
IS8HB06	Blue	05391541270454	05391541271024
IS8HB06W	White	05391541272854	05391541272861
IS8HB06G	Green	05391541272878	05391541272885
IS8HB06Y	Yellow	05391541272892	05391541272908



ISOL8 theatre headwear products are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 16/05/25, Version 004