

# ISOL/8<sup>®</sup>

## SPECIFICATION SHEET

*Designed by Healthcare Professionals,  
For Healthcare Professionals*



### **Basics Theatre Hijab**

File no.: SPEC-IS8HB07

#### **ISOL8 HEALTHCARE LTD**

**Head office:**

Unit 4 Joyce House,  
Barrack Square, Ballincollig,  
Co Cork, P31 HW35,  
Ireland

**UK office:**

Unit 5 Catkin Way, St Helen  
Auckland, Bishop Auckland,  
Co Durham, DL14 9TF,  
United Kingdom

#### **CONTACT**

info@isol8health.com  
IE: +353 (0) 66 71 21659  
GB: +44 (0) 191 816 1966  
www.isol8health.com

# Complete coverage of the head, neck and shoulders

Our range of Basics headwear are all made from a soft, breathable spunbond polypropylene fabric which offers a great combination of both comfort and affordability. Our Basics range has been designed to meet the basic needs of your operating department while offering great value for your money.

This single use theatre hijab offers our greatest level of cover by completely enshrouding the head, neck and shoulders. An elasticated front tucks under the wearers chin exposing only the face. It is an ideal product for procedures where extra cover is required or for religious reasons.



**Basics**

## Dimensions

Crown (Length x Width)	22 x 21 cm $\pm$ 1cm
Depth (front)	5cm $\pm$ 1cm
Depth (rear)	54cm $\pm$ 3cm

# Product Specification

Product Name	Basics Theatre Hijab
Manufacturer Product Codes (MPN)	IS8HB07
NPC	BWF85018
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Single use theatre hijab offering full coverage across the head, neck and shoulders
Sizes (see dimensions below)	One size fits all
Fabric	Spunbond Polypropylene 25gsm
Colour	Blue
Sterility	Non-sterile
Ties	Elasticated front exposing the face only
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner box)	40 hijabs per dispenser box 330x140x185 mm
Product packaging (outer carton)	6 boxes per carton (240 hijabs)
Recycling Information	Outer carton and inner box are 100% recyclable (cardboard) Hats are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541271062
GTIN code outer carton	05391541271116





ISOL8 theatre headwear products are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD  
Unit 4 Joyce House, Barrack Square, Ballincollig,  
Co Cork, P31 HW35. Company No. 687076

[info@isol8health.com](mailto:info@isol8health.com)  
IE: +353 (0) 66 71 21659  
GB: +44 (0) 191 816 1966  
[www.isol8health.com](http://www.isol8health.com)

# ISOL/8<sup>®</sup>

**Designed by Healthcare Professionals,  
For Healthcare Professionals**

Date of issue: 16/05/25, Version 004