

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Short-sleeved Patient Examination Gown

File no.: SPEC-IS8PG

ISOL8 HEALTHCARE LTD

Head office:

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland

UK office:


Unit 5 Catkin Way, St Helen
Auckland, Bishop Auckland,
Co Durham, DL14 9TF,
United Kingdom

CONTACT


info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

Comfortable and high level of dignity

Our single-use short-sleeved patient examination gown is a soft breathable gown which is opaque in colour to offer the wearer a high level of dignity. It closes fully at the rear by way of 2 ties and allows easy examination of the patient by healthcare professionals during patient encounters and surgical or diagnostic procedures. The SMS fabric offers protection against light splashes and mechanical movements. This gown comes in size large which is built to fit most body shapes and sizes and also comes in size XXL for large framed persons.



Product Photos



Product Specification

Product Name	Short-sleeved Patient Examination Gown
Manufacturer Product Codes (MPN)	IS8PG01 (L), IS8PG02 (XXL)
Manufacturer	ISOL8 Healthcare Ltd.
Gown type	Opaque patient examination gown – Single use
Sizes (see size guide below for dimensions)	L (Suitable for most sizes) XXL (Large framed persons)
Fabric	Extra soft SMS 48gsm (100% polypropylene)
Gown cut	Short sleeved
Neck closure	Ties
Colour	Navy blue
Sterility	Non-sterile
Ties	Rear ties
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner pack)	10 gowns per (labelled) polythene bag
Product packaging (outer carton)	10 x packs per carton (100 gowns) Outer carton: 600x390x360 mm (L), 600x390x380 mm (XXL)
Recycling Information	Outer carton is 100% recyclable (cardboard), Inner pack is 100% recyclable (PE), Gowns are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485

Size Guide

Size	Length (cm)	Hem circumference (cm)
L	105	150
XXL	105	165

Product Codes

MPN	NPC	Size	Quantity/carton	GTIN (gown)	GTIN (carton)
IS8PG01	BWK85007	L	100	05391541270515	05391541270904
IS8PG02	BWK85008	XXL	100	05391541270522	05391541270911



ISOL8 patient gowns are CE certified class I medical devices intended for hospital patients. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2025 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 14/08/25, Version 009