

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



FLUIDSAFE Child-friendly Procedure Mask

File no.: SPEC-IS8RM04P

ISOL8 HEALTHCARE LTD

Head office:

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland

UK office:

Unit 5 Catkin Way, St Helen
Auckland, Bishop Auckland,
Co Durham, DL14 9TF,
United Kingdom

CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

Fluid resistant and extra soft

Our **FLUIDSAFE** standard procedure mask is manufactured from high quality fabric with a super soft silk-like inner material providing excellent comfort to the wearer. The fun and exciting pattern of underwater creatures is loved by young children and is perfect for relieving their anxiety in paediatric theatres and departments. Our type IIR masks provide the highest barrier performance and their fluid resistant properties offer protection against splashing during procedures.

Our procedure masks are most suited for use in clinical environments for patient encounters, examinations and minor procedures. Elastic earloops are easy to apply and provide a comfortable fit.



Packaging Information

Product packaging (inner pack)	50 masks per dispenser box Dispenser box: 190x102x85 mm
Product packaging (outer carton)	40 inner boxes per carton (2000 masks), Outer carton: 430x400x440 mm, Gross weight: 13Kg
Recycling Information	Inner box and outer carton are 100% recyclable (cardboard), Mask is non-recyclable

Product Specification

Product Name	FLUIDSAFE Child-friendly Procedure Mask
Manufacturer Product Codes (MPN)	IS8RM04P
NPC	BWG85083
Manufacturer	ISOL8 Healthcare Ltd.
Mask type	Fluid resistant child-friendly adult mask for use by paediatric healthcare professionals
Compliance with standard	EN14683:2019 Type IIR
Mask size (Length x Width)	175x95 mm
Mask fit	Ties
Colour	White with child-friendly animal pattern
Sterility	Non-sterile
Mask fabric (3-layers)	Outer: 35gsm Spunbond Polypropylene, Middle: 30gsm Melt-blown Polypropylene, Inner: 18gsm extra soft ES silk-like fibre
Earloops	Spandex polyurethane coated with nylon thread 160±5mm x 3mm
Nose bar material	Full length aluminium nose piece
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541270348
GTIN code outer carton	05391541270874



ISOL8 masks are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**