

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



SMS Medical Gown

File no.: SPEC-IS8SG

ISOL8 HEALTHCARE LTD

Unit 4 Joyce House,
Barrack Square, Ballincollig,
Co Cork, P31 HW35,
Ireland


CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

Tri-layer SMS for high level of protection and comfort

Our single-use SMS medical gown is a soft, breathable gown designed and tested to be resistant against bodily fluids, microorganisms and contaminated solids. This tri-layer SMS gown ensures a high level of protection while keeping the wearer comfortable during patient encounters and non-sterile procedures.

This gown comes in size large which is built to fit most body shapes and sizes and also comes in size XXL for large framed persons.



Product Photos



Product Specification

Product Name	SMS Medical Gown
Manufacturer Product Codes (MPN)	IS8SG01 (L), IS8SG02 (XXL)
Manufacturer	ISOL8 Healthcare Ltd.
Gown type	Non-sterile medical gown – Single use
Compliance with standard	EN13795-1:2019 - Standard Performance
Sizes (see size chart below for dimensions)	L (Suitable for most sizes) XXL (Large framed persons)
Fabric	SMS 40gsm (100% polypropylene)
Gown cut	Standard cut
Sleeve Seam	Double overlock stitching
Neck closure	Hook and loop
Cuffs	Soft knitted cuffs
Colour	Blue
Sterility	Non-sterile
Ties	Ultrasonically welded waistbelt (tie at rear)
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	5 years
Product packaging (inner pack)	10 gowns per (labelled) polythene bag
Product packaging (outer carton)	10 x packs per carton (100 gowns) Outer carton: 600x390x400 mm (L), 600x390x420 mm (XXL)
Recycling Information	Outer carton is 100% recyclable (cardboard), Inner pack is 100% recyclable (PE), Gowns are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485

Size Guide

Size	Length (cm)	Hem circumference (cm)
L	120	155
XXL	130	165

Product Codes

MPN	NPC	Size	Quantity/carton	GTIN (gown)	GTIN (carton)
IS8SG01	BWK85009	L	100	05391541270492	05391541270881
IS8SG02	BWK85010	XXL	100	05391541270508	05391541270898



ISOL8 medical gowns are CE certified class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2026 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 21/01/26, Version 006